

# Acquisition Protocol for Cytex<sup>®</sup> cFluor<sup>®</sup> Human B Cell Monitoring Kit

## Contents

Introduction .....	1
Preparing SpectroFlo <sup>®</sup> Software .....	2
Add fluorochromes to the library .....	2
Import the experiment template .....	2
Setting up the Instrument .....	2
Acquiring Controls and Samples .....	4
Acquire controls in the Reference Group .....	4
Unmix reference controls .....	5
Acquire and analyze multicolor samples .....	6
Appendix A: Acquisition Setup .....	10
Appendix B: Reusing Single Color Controls .....	11
Appendix C: Single Color Control Gating and Signatures for 3-Laser (V-B-R) Cytex <sup>®</sup> Cytometers .....	12
Appendix D: Example of Similarity <sup>™</sup> Matrix for Cytex <sup>®</sup> cFluor <sup>®</sup> Human B Cell Monitoring Kit .....	21
Appendix E: Example of Gating Multicolor Samples .....	23

## Introduction

This acquisition protocol provides step-by-step instructions to set up your Cytex<sup>®</sup> Northern Lights<sup>®</sup> system (2-laser B-R configuration or higher) or Aurora<sup>®</sup> system for data acquisition of the Cytex<sup>®</sup> cFluor<sup>®</sup> Human B Cell Monitoring Kit. This protocol provides instructions on 1) preparing SpectroFlo<sup>®</sup> software, 2) setting up the instrument, and 3) acquiring controls and samples.

This kit contains 25 tests.

**NOTE:** For suggestions on how to prepare human peripheral mononuclear cells and whole blood, see **Sample Preparation (PBMCs) Guidelines for Cytex<sup>®</sup> cFluor<sup>®</sup> Human B Cell Monitoring Kit** and **Sample Preparation (Whole Blood) Guidelines for Cytex<sup>®</sup> cFluor<sup>®</sup> Human B Cell Monitoring Kit**.


## Preparing SpectroFlo® Software

### Add fluorochromes to the library

1. Add the kit fluorochromes into the **Library**. Select **Fluorescent Tags** in the **Library** and import the **Cytek® 13C B cell Assay Tags.csv** file.

**NOTE:** If any fluorochromes were previously entered into SpectroFlo® Library, a warning message will appear. Click **Replace** to overwrite the information in the library with the new fluorescent tag information.

Importing Fluorescent Tags:




Name	Emission Wavelength	Laser Excitation	Display Name	Synonyms
cFluor B515	515	488	cFluor B515	cFluorB515
cFluor B548	548	488	cFluor B548	cFluorB548
cFluor BYG575	575	488	cFluor BYG575	cFluorBYG575
cFluor BYG667	667	488	cFluor BYG667	cFluorBYG667
cFluor BYG710	710	488	cFluor BYG710	cFluorBYG710
cFluor BYG781	781	488	cFluor BYG781	cFluorBYG781
cFluor R659	659	635	cFluor R659	cFluorR659
cFluor R668	668	635	cFluor R668	cFluorR668
cFluor R685	685	635	cFluor R685	cFluorR685
cFluor R720	720	635	cFluor R720	cFluorR720
cFluor R780	780	635	cFluor R780	cFluorR780
cFluor B532	532	488	cFluor B532	cFluorB532
cFluor BYG610	610	488	cFluor BYG610	cFluorBYG610
ViaDye Red	740	635	ViaDye Red	

### Import the experiment template

1. The **Cytek® B Panel Assay Templates** include a reference group with predefined stopping criteria, assigned marker names, as well as recommended acquisition and analysis worksheets. The Template for PBMCs is setup for ViaDye Red, while the Blood Template is not. The use of viability dyes is not needed with fresh blood.
2. Add Templates into the **Library**. Select **Experiment Templates** and import the templates.

Importing Experiment Templates:




Name	Created By	Date Created	Description
B cell Panel QC Template	Admin	October 28, 2022 - 08:57 AM	
Cytek B Panel Assay Template_Blood	Admin	October 05, 2022 - 10:50 AM	
Cytek B Panel Assay Template_PBMC	Admin	October 05, 2022 - 10:09 AM	
Default	Admin	July 13, 2021 - 09:10 AM	Default Experiment

## Setting up the Instrument

1. Follow the instructions for instrument setup and **Performing Daily QC** as outlined in the User's Guide.

- From the **Acquisition** module, create a new experiment by clicking on **Template**, then choose either “Cytek® B Panel Assay Template\_PBMC” or “Cytek® B Panel Assay Template\_Blood”

Opening Experiment Template:



- Either beads or cells can be used for single color controls except cFluor® R668 IgG for which only beads should be used. See below (Table 1) for reference control sample type recommendations for each marker.

**Table 1.** Reference Control Type Recommendations


Laser	Target	Fluorochrome	Recommended Control Type
Blue	IgM	cFluor® B515	Cells or Beads
	CD4	cFluor® B532	Cells or Beads
	CD15	cFluor® B548	Cells or Beads
	CD38	cFluor® BYG575	Cells or Beads
	CD27	cFluor® BYG610	Cells or Beads
	CD14	cFluor® BYG667	Cells or Beads
	IgD	cFluor® BYG710	Cells or Beads
	CD19	cFluor® BYG781	Cells or Beads
Red	CD3	cFluor® R659	Cells or Beads
	IgG	cFluor® R668	Beads*
	CD8	cFluor® R685	Cells or Beads
	CD20	cFluor® R720	Cells or Beads
	CD45	cFluor® R780	Cells or Beads

\*If cells must be used for the single-color control of cFluor R668 anti-IgG, since the IgG+ B cell population is small, it is recommended to use either cFluor R668 conjugated anti-human CD3 or CD4, or other antibodies for highly expressed markers.

- The Experiment Template is setup for using cells as references. If using beads, go to **Edit**, select the **Groups** tab, and click on **Reference Group**. In **Edit Reference Group**, for the Fluorescent Tags that will be using beads, change the **Control Type** to “Beads.” Be sure to save all changes.

- Each sample tube is set to acquire a certain number of cells (see Appendix A). This can be changed by clicking **Edit**, then changing the Stopping Gate and Stopping Criteria under **Acquisition**. If using beads as references, the stopping criteria should be changed to collect 5,000 beads for all reference tubes. After making any changes to the acquisition criteria be sure to save all changes.
- Add tubes and groups as needed for acquiring all samples. To preserve the predefined acquisition conditions, duplicate the existing tubes or groups.

Changing Reference Controls from Cells to Beads in Experiment Template:



The screenshot shows the Cytek software interface for editing an experiment template. The left sidebar displays a tree view of the experiment setup, including a 'Reference Group' and a 'MultiColor' group. The top toolbar contains buttons for 'Save', 'Save As', 'Edit', and 'Reference Group'. The main window is divided into two tabs: 'Edit Experiment' and 'Edit Reference Group'. The 'Edit Reference Group' tab is active, showing a table of fluorescent tags and their control types. The 'Control Type' column is highlighted, and 'Beads' is selected for the 'Control Type' of the 'cFluor B515' tag. The 'Save' button is highlighted in the bottom right corner.


Fluorescent Tag	Control Type	Label	Negative Control
cFluor B515	Beads	IgM	
cFluor B532	Beads	CD4	
cFluor B548	Beads	CD15	
cFluor BYG575	Beads	CD38	
cFluor BYG610	Beads	CD27	
cFluor BYG667	Beads	CD14	
cFluor BYG710	Beads	IgD	
cFluor BYG781	Beads	CD19	
cFluor R659	Beads	CD3	
cFluor R668	Beads	CD3 (IgG)	
cFluor R685	Beads	CD8	
cFluor R720	Beads	CD20	
cFluor R780	Beads	CD45	

## Acquiring Controls and Samples

### Acquire controls in the Reference Group

- Preview unstained cell control at low flow rate to minimize wasted sample volume. Starting from the default CytekAssaySetting, optimize the FSC and SSC gains, as well as the threshold to fully visualize the cells of interest (see Figure 1).

**NOTE:** Instrument settings can be saved as "Cytek® B cell Kit" for future use by clicking the **Save As** button in the Instrument Control window. The gains for all fluorescent parameters are set up with CytekAssaySetting in the instrument and only FCS, SSC and threshold need to be optimized for specific sample types.




**Figure 1:** Example cell resolution of human blood samples after FSC and SSC adjustment with threshold set at FSC 350,000.

2. Acquire unstained and single-color controls in the Reference Control group. Click **Start** to preview for 5 to 10 seconds until the event rate stabilizes, then click **Record** to record each sample.

### Unmix reference controls

1. Once all controls have been acquired, Click **Unmix**.  
**NOTE:** Refer to Appendix B for additional workflows to reuse the reference controls.
2. Under **Select Controls** tab in the Unmix Experiment wizard, ensure:
  - 1) Under **Unstained Controls** that Unstained (Cells) is selected as the Reference Group
  - 2) Under **Stained Controls** that the Unstained column is blank, which will enable selecting the negative population in the **Identify Positive/Negative Populations** section
  - 3) **Autofluorescence as a Fluorescence Tag** is selected for 5 laser instruments and unselected for 2-4 laser instruments. Click **Next**.

### Selecting Controls for Unmixing:



**Unmix Experiment**

Select Controls   Identify Positive/Negative Populations   QC Controls

**UNSTAINED CONTROLS**

☐ Use Control from Library

☒ Use Control from Experiment   Reference Group - Unstained (Cells)

Name	Control Type
Reference Group - Unstained (Cells)	Cells

**STAINED CONTROLS**

<input type="checkbox"/> From Library	Fluorescent Tag	Control	Unstained	Generic
<input type="checkbox"/>	cFluor B515	IgM cFluor B515 (Cells)	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
<input type="checkbox"/>	cFluor B532	CD4 cFluor B532 (Cells)	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
<input type="checkbox"/>	cFluor B548	CD15 cFluor B548 (Cells)	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
<input type="checkbox"/>	cFluor BYG575	CD38 cFluor BYG575 (Cells)	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
<input type="checkbox"/>	cFluor BYG610	CD27 cFluor BYG610 (Cells)	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
<input type="checkbox"/>	cFluor BYG667	CD14 cFluor BYG667 (Cells)	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
<input type="checkbox"/>	cFluor BYG710	IgD cFluor BYG710 (Cells)	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
<input type="checkbox"/>	cFluor BYG781	CD19 cFluor BYG781 (Cells)	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>

☐ Autofluorescence as a Fluorescent Tag   Only apply for 5 laser instruments

- Under **Identify Positive/Negative Populations** tab, ensure for all single color controls that:

- 1) Scatter plot is gated on the appropriate populations,
- 2) Black bar is on the peak channel,
- 3) Signature of each fluorochrome matches the expected spectrum, and
- 4) Negative and positive gates in the histogram are correctly positioned.


**NOTE:** Refer to Appendix C for the correct gate positioning, expected spectra and peak channels of each fluorochrome, and the positioning of negative and positive gates in the histograms.

3. Click **Next**. Under the **QC Controls** tab, click on **Similarity Matrix** to confirm all controls were appropriately stained. Click on **View Similarity Index** to compare the expected complexity index value found in Appendix D.
4. Click **Live Unmixing**.

### Acquire and analyze multicolor samples

1. Acquire multicolor samples. See Appendix E for example of gating populations.
2. Compensation should not be needed when using cells as references for unmixing on a 2 or 3 laser system. There may be minor unmix errors using beads and in the BYG dyes when using a 5-laser system, which usually does not affect the identification of cell populations of interest. To adjust compensation, open the **B Panel\_1xN plot** worksheet.
3. Right click on the multicolor tube needing compensation adjustments. From the drop-down menu, click **Edit Properties**.
4. Click on **Enable Compensation** in the pop-up wizard. Leave the wizard open and move it aside.

## Opening Spillover Matrix to Adjust Compensation if Needed:



The screenshot shows the 'Tube/Well Properties' dialog box with the 'Spillover' tab selected. The 'Enable Compensation' checkbox is checked. The spillover matrix table is displayed below.


INTO → FROM →	cFluor B515	cFluor B532	cFluor B548	cFluor BYG575	cFluor BYG610	cFluor BYG667	cFluor BYG710	cFluor BYG781	cFluor R659	cFluor R668	cFluor R685	cFluor R720
cFluor B515	100	0	0	0	0	0	0	0	0	0	0	0
cFluor B532	0	100	0	0	0	0	0	0	0	0	0	0
cFluor B548	0	0	100	0	0	0	0	0	0	0	0	0
cFluor BYG575	0	0	0	100	0	0	0	0	0	0	0	0
cFluor BYG610	0	0	0	0	100	0	0	0	0	0	0	0
cFluor BYG667	0	0	0	0	0	100	0	0	0	0	0	0
cFluor BYG710	0	0	0	0	0	0	100	0	0	0	0	0
cFluor BYG781	0	0	0	0	0	0	0	100	0	0	0	0
cFluor R659	0	0	0	0	0	0	0	0	100	0	0	0
cFluor R668	0	0	0	0	0	0	0	0	0	100	0	0
cFluor R685	0	0	0	0	0	0	0	0	0	0	100	0
cFluor R720	0	0	0	0	0	0	0	0	0	0	0	100

### 5. To check and adjust for unmixing error:

- 1) Select all permutation plots,
- 2) Right click and select **Properties** from the drop-down menu
- 3) From **Plot Properties**, select the first fluorochrome under **X Axis Parameter**, and
- 4) Check all permutation plots against the first fluorochrome for any unmixing errors and adjust spillover (compensation) as needed. Do this for all fluorochromes by selecting each fluor one by one under X Axis Parameter.




## Using 1XN Plots to Assess Unmixing Errors:



- To adjust unmixing error of a plot, click the **Adjust Spillover** icon in the ribbon menu. Click and drag upward or downward on the plot to make the adjustment (see Figure 2). Alternatively, the spillover values can be typed directly into the Spillover Matrix by double clicking on the cell you want to change the value for.

## Adjusting Spillover to Correct Unmix Error:




**Figure 2:** Example of unmixing error adjustment for a blood sample unmixed with beads.



7. Click **Save** and close out of the **Tube/Well Properties** wizard. If similar unmix errors are observed among the multicolor samples, the adjusted Spillover matrix can be applied from one Tube to other Tubes by copying and pasting the Spillover into those Tubes. To do this, right click on Tubes and use the drop down menu to **Copy Spillover** and **Paste Spillover**.
8. Repeat the unmixing error adjustment for all multicolor samples as needed.
9. The manually adjusted compensation matrix can be saved to the **Library** for future use if using SpectroFlo® version 3.0 or higher. To do this, open the **Tube/Well Properties** and in the **Spillover** section click the **Save to Library** icon in the ribbon menu directly above the spillover matrix.
10. For convenience of analyzing samples quickly, we provide an Analysis worksheet that has data plots with gating for populations of interest and a corresponding statistics table with cell counts and percentages.
11. To export an experiment:
  - 1) Save and close the experiment,
  - 2) Click **My Experiments**,
  - 3) Select the experiment and click **Export**, and
  - 4) Choose a directory and click **Export**.

#### Exporting Experiment:



## Appendix A: Acquisition Setup

Example of Acquisition Setup for PBMCs using cells as References. For Acquiring Blood, the setup does not include ViaDye Red and for acquiring multicolor samples the Events to Record is 1 million Cells. For using Beads as References, the Events To Record can be set for acquiring 5,000 Beads for each control.

Edit Experiment

Fluorescent Tags
Groups
Markers
Keywords
Acquisition

☐ Tube/Well Specific User Setting
Experiment User Setting:
CytekAssaySetting (Cytek)

Name	Worksheet	Stopping Gate	Storage Gate	Events To Record	Stopping Volume (ul)	Stopping Criteria	Stopping Time (sec)
▼ B cell Kit_PBMC				1 - 10,000,000	3,000	<input type="checkbox"/> Count & Volume	10,000
▼ Reference Group	B_Panel_References (Raw)	▼ Lymph or Beads	▼ All Events	▼ 1 - 10,000,000	3,000	<input type="checkbox"/> Count & Volume	10,000
Unstained (Cells)	B_Panel_References (Raw)	▼ Lymph or Beads	▼ All Events	▼ 5,000	3,000	<input type="checkbox"/> Count & Volume	10,000
IgM cFluor B515 (Cells)	B_Panel_References (Raw)	▼ Lymph or Beads	▼ All Events	▼ 10,000	3,000	<input type="checkbox"/> Count & Volume	10,000
CD4 cFluor B532 (Cells)	B_Panel_References (Raw)	▼ Lymph or Beads	▼ All Events	▼ 5,000	3,000	<input type="checkbox"/> Count & Volume	10,000
CD4 (CD15) cFluor B548 (Cells)	B_Panel_References (Raw)	▼ Lymph or Beads	▼ All Events	▼ 10,000	3,000	<input type="checkbox"/> Count & Volume	10,000
CD38 cFluor BYG575 (Cells)	B_Panel_References (Raw)	▼ Lymph or Beads	▼ All Events	▼ 10,000	3,000	<input type="checkbox"/> Count & Volume	10,000
CD27 cFluor BYG610 (Cells)	B_Panel_References (Raw)	▼ Lymph or Beads	▼ All Events	▼ 5,000	3,000	<input type="checkbox"/> Count & Volume	10,000
CD14 cFluor BYG667 (Cells)	B_Panel_References (Raw)	▼ Lymph or Beads	▼ All Events	▼ 10,000	3,000	<input type="checkbox"/> Count & Volume	10,000
IgD cFluor BYG710 (Cells)	B_Panel_References (Raw)	▼ Lymph or Beads	▼ All Events	▼ 30,000	3,000	<input type="checkbox"/> Count & Volume	10,000
CD19 cFluor BYG781 (Cells)	B_Panel_References (Raw)	▼ Lymph or Beads	▼ All Events	▼ 10,000	3,000	<input type="checkbox"/> Count & Volume	10,000
CD3 cFluor R659 (Cells)	B_Panel_References (Raw)	▼ Lymph or Beads	▼ All Events	▼ 5,000	3,000	<input type="checkbox"/> Count & Volume	10,000
CD3 (IgG) cFluor R668 (Cells)	B_Panel_References (Raw)	▼ Lymph or Beads	▼ All Events	▼ 5,000	3,000	<input type="checkbox"/> Count & Volume	10,000
CD8 cFluor R685 (Cells)	B_Panel_References (Raw)	▼ Lymph or Beads	▼ All Events	▼ 10,000	3,000	<input type="checkbox"/> Count & Volume	10,000
CD20 cFluor R720 (Cells)	B_Panel_References (Raw)	▼ Lymph or Beads	▼ All Events	▼ 10,000	3,000	<input type="checkbox"/> Count & Volume	10,000
ViaDye Red (Cells)	B_Panel_References (Raw)	▼ Lymph or Beads	▼ All Events	▼ 5,000	3,000	<input type="checkbox"/> Count & Volume	10,000
CD45 cFluor R780 (Cells)	B_Panel_References (Raw)	▼ Lymph or Beads	▼ All Events	▼ 5,000	3,000	<input type="checkbox"/> Count & Volume	10,000
▼ MultiColor	B_Panel_PBMC_Acquire (Unmixed)	▼ Cells	▼ All Events	▼ 500,000	3,000	<input type="checkbox"/> Count & Volume	10,000
Tube1	B_Panel_PBMC_Acquire (Unmixed)	▼ Cells	▼ All Events	▼ 500,000	3,000	<input type="checkbox"/> Count & Volume	10,000
Tube2	B_Panel_PBMC_Acquire (Unmixed)	▼ Cells	▼ All Events	▼ 500,000	3,000	<input type="checkbox"/> Count & Volume	10,000

Previous
Cancel
Save and Open


## Appendix B: Reusing Single Color Controls

To reuse reference controls in future experiments, follow the instructions below.

**NOTE:** For best results, maintain the instrument properly, perform QC daily calibration and use the same reagent lot. The samples need to be collected using the same instrument as the reference controls.

### Reusing single color controls from a previous experiment

1. From the **Acquisition** module, open **My Experiments**. Right click the saved experiment and select **Duplicate with Reference Group Data**. This will duplicate the experiment with the reference controls and multicolor samples.




2. Open the newly created Experiment
3. Under Multicolor Group, add tubes for the new multicolor samples as needed
4. Click **Edit** to open Edit Experiment wizard
5. Under **Acquisition**, check to make sure the settings have the correct Worksheets, Stopping Gates, and Events to Record
6. Click **Save and Open**
7. Preview a sample to set FSC and SSC gains and FSC threshold
8. Acquire multicolor samples
9. Adjust the gates in the analysis worksheet as needed.


## Appendix C: Single Color Control Gating and Signatures for 3-Laser (V-B-R)


### Cytek® Cytometers

Example of gating and spectrum signatures for single color controls using blood.














Example of gating and spectrum signatures for single color controls using PBMCs.












Example of gating and spectrum signatures for Cytex® FSP™ CompBeads.










## Appendix D: Example of Similarity™ Matrix for Cytex® cFluor® Human B Cell Monitoring Kit

Expected Similarity™ matrix and Complexity™ index generated on a 3-laser (V-B-R) Cytex® Aurora®. Generating similar values is a good indication that signatures of your single color controls match those generated by Cytex®.


Similarity matrix and Complexity index generated from blood sample:






## Appendix E: Example of Gating Multicolor Samples

Gating of Populations for Blood acquired on 3 laser Aurora®



Gating of Populations for PBMC acquired on 3 laser Aurora®





**For Research Use Only. Not intended for use in diagnostic procedures.**

Cytek® FSP™ CompBeads are developed and manufactured by Slingshot Biosciences, Inc.

cFluor® B515, cFluor® B532, cFluor® B548, cFluor® R668 cFluor® R685 and cFluor® R720 are equivalent to CF® 488A, CF® 503, CF® 514, CF® 647, CF® 660C and CF® 700 respectively, manufactured and provided by Biotium, Inc. under an Agreement between Biotium and Cytek (LICENSEE). The manufacture, use, sale, offer for sale, or import of the product is covered by one or more of the patents or pending applications owned or licensed by Biotium. The purchase of this product includes a limited, non-transferable immunity from suit under the foregoing patent claims for using only this amount of product for the purchaser's own internal research. No right under any other patent claim, no right to perform any patented method, and no right to perform commercial services of any kind, including without limitation reporting the results of purchaser's activities for a fee or other commercial consideration, is conveyed expressly, by implication, or by estoppel.

"Cytek", "Aurora", "Northern Lights", "FSP", "cFluor" and "ViaDye" are trademarks or registered trademarks of Cytek Biosciences, Inc. All other service marks, trademarks and tradenames appearing herein are the property of their respective owners.