

## APC/Fire™ 810 anti-human CD11b Antibody

|                          |  |
|--------------------------|--|
| <b>Catalog# / Size</b>   | 301371 / 25 tests<br>301372 / 100 tests  |
| <b>Clone</b>             | ICRF44   |
| <b>Regulatory Status</b> | RUO  |
| <b>Workshop</b>          | IV M047  |
| <b>Other Names</b>       | Integrin $\alpha_M$ chain, C3biR, CR3, Mac-1, Mo1, ITGAM   |
| <b>Isotype</b>           | Mouse IgG1, $\kappa$   |
| <b>Description</b>       | CD11b is a 165-170 kD type I transmembrane glycoprotein also known as $\alpha_M$ integrin, Mac-1, CR3, and C3biR. CD11b non-covalently associates with integrin $\beta_2$ (CD18) and is expressed on granulocytes, monocytes/macrophages, dendritic cells, NK cells, and subsets of T and B cells. CD11b/CD18 is critical for the transendothelial migration of monocytes and neutrophils. It is also involved in granulocyte adhesion, phagocytosis, and neutrophil activation. CD11b/CD18 interacts with ICAM-1 (CD54), ICAM-2 (CD102), ICAM-4, CD14, CD23, heparin, iC3b, fibrinogen, and factor X. |

### Product Details

---

|                               |   |
|-------------------------------|---|
| <b>Verified Reactivity</b>    | Human, Cynomolgus, Rhesus   |
| <b>Reported Reactivity</b>    | African Green, Baboon, Chimpanzee, Common Marmoset, Pig   |
| <b>Antibody Type</b>          | Monoclonal  |
| <b>Host Species</b>           | Mouse   |
| <b>Formulation</b>            | Phosphate-buffered solution, pH 7.2, containing 0.09% sodium azide and BSA (origin USA)   |
| <b>Preparation</b>            | The antibody was purified by affinity chromatography and conjugated with APC/Fire™ 810 under optimal conditions.  |
| <b>Concentration</b>          | Lot-specific (to obtain lot-specific concentration and expiration, please enter the lot number in our <a href="#">Certificate of Analysis</a> online tool.)   |
| <b>Storage &amp; Handling</b> | The antibody solution should be stored undiluted between 2°C and 8°C, and protected from prolonged exposure to light. <b>Do not freeze.</b>   |
| <b>Application</b>            | <a href="#">FC - Quality tested</a>   |
| <b>Recommended Usage</b>      | <p>Each lot of this antibody is quality control tested by <a href="#">immunofluorescent staining with flow cytometric analysis</a>. For flow cytometric staining, the suggested use of this reagent is 5 <math>\mu</math>L per million cells in 100 <math>\mu</math>L staining volume or 5 <math>\mu</math>L per 100 <math>\mu</math>L of whole blood. It is recommended that the reagent be titrated for optimal performance for each application.</p> <p>* APC/Fire™ 810 has a maximum excitation of 650 nm and a maximum emission of 810 nm.</p> <p>Excessive exposure to light, and commonly used fixation, permeabilization buffers can affect APC/Fire™ 810 fluorescence signal intensity and spread. Please keep conjugates protected from light exposure. For more information and representative data, visit our <a href="#">Fire Dyes</a> page.</p> |
| <b>Excitation Laser</b>       | Red Laser (633 nm)  |
| <b>Application Notes</b>      | The ICRF44 antibody inhibits heterotypic adhesion of granulocytes in response to fMLP. Additional reported applications (for the relevant formats) include: immunohistochemical staining of acetone-fixed frozen tissue sections, immunofluorescence microscopy <sup>5</sup> , stimulation of monocytes <sup>3</sup> , blocking of heterotypic PMN aggregation <sup>8</sup> , and blocking of granulocyte activation <sup>12</sup> . This clone was tested in-house and does not work on formalin fixed paraffin-embedded (FFPE) tissue.  |

The Ultra-LEAF™ purified antibody (Endotoxin < 0.01 EU/μg, Azide-Free, 0.2 μm filtered) is recommended for functional assays (Cat. Nos. 301361 & 301362).

#### Application References

(PubMed link indicates BioLegend citation)

1. Knapp W. 1989. Leucocyte Typing IV. Oxford University Press New York.
2. Barclay N, *et al.* 1997. The Leucocyte Antigen Facts Book. Academic Press Inc. San Diego.
3. Rezzonico R, *et al.* 2001. *Blood* 97:2932. (Stim)
4. Marsik C, *et al.* 2003. *Shock* 20:493. (FC)
5. David A, *et al.* 2003. *J. Leukoc. Biol.* 74:551. (IF)
6. Charles N, *et al.* 2010. *Nat. Med.* 16:701. (FC) [PubMed](#)
7. Thurlow LR, *et al.* 2010. *Infect. Immun.* 128:1128. (FC) [PubMed](#)
8. Jadhav S, *et al.* 2001. *J. Immunol.* 167:5986. (Block)
9. Yoshino N, *et al.* 2000. *Exp. Anim. (Tokyo)* 49:97. (FC)
10. Sestak K, *et al.* 2007. *Vet. Immunol. Immunopathol.* 119:21. (FC)
11. Wen T, *et al.* 2014. *J Immunol.* 192:5481. (FC) [PubMed](#)
12. Sprong T, *et al.* 2003. *Blood* 102:3702. (Block)

[See More](#)

**RRID** AB\_3683301 (BioLegend Cat. No. 301371)  
AB\_3683301 (BioLegend Cat. No. 301372)

## Antigen Details

---

|                           |  |
|---------------------------|--|
| <b>Structure</b>          | Integrin, type I transmembrane glycoprotein, associates with integrin $\beta_2$ (CD18), 165-170 kD                         |
| <b>Distribution</b>       | Granulocytes, monocytes/macrophages, dendritic cells, NK cells, subset of T cells, subset of B cells                       |
| <b>Function</b>           | Adhesion, phagocytosis, chemotaxis, neutrophil activation  |
| <b>Ligand/Receptor</b>    | ICAM-1(CD54), ICAM-2 (CD102), ICAM-4, CD14, CD23, heparin, iC3b, fibrinogen, factor X                                      |
| <b>Cell Type</b>          | B cells, Dendritic cells, Granulocytes, Macrophages, Monocytes, Neutrophils, NK cells, T cells, Tregs                      |
| <b>Biology Area</b>       | Cell Adhesion, Cell Biology, Costimulatory Molecules, Immunology, Innate Immunity, Neuroscience, Neuroscience Cell Markers |
| <b>Molecular Family</b>   | Adhesion Molecules, CD Molecules   |
| <b>Antigen References</b> | 1. Stewart M, <i>et al.</i> 1995. <i>Curr Opin Cell Biol.</i> 7:690.   |
| <b>Gene ID</b>            | <a href="#">3684</a>   |

## Related Protocols

---

- [Cell Surface Flow Cytometry Staining Protocol](#)


## Other Formats

---

APC anti-human CD11b, Biotin anti-human CD11b, PE anti-human CD11b, PE/Cyanine5 anti-human CD11b, Purified anti-human CD11b, Pacific Blue™ anti-human CD11b, Alexa Fluor® 488 anti-human CD11b, Alexa Fluor® 647 anti-human CD11b, PE/Cyanine7 anti-human CD11b, PerCP/Cyanine5.5 anti-human CD11b, Brilliant Violet 421™ anti-human CD11b, Brilliant Violet 570™ anti-human CD11b, FITC anti-human CD11b, Brilliant Violet 605™ anti-human CD11b, Brilliant Violet 510™ anti-human CD11b, Brilliant Violet 650™ anti-human CD11b, Purified anti-human CD11b (Maxpar® Ready), Alexa Fluor® 594 anti-human CD11b, APC/Cyanine7 anti-human CD11b, Brilliant Violet 711™ anti-human CD11b, Brilliant Violet 785™ anti-human CD11b, PE/Dazzle™ 594 anti-human CD11b, APC/Fire™ 750 anti-human CD11b, TotalSeq™-A0161 anti-human CD11b, Alexa Fluor® 700 anti-human CD11b, TotalSeq™-B0161 anti-human CD11b, TotalSeq™-C0161 anti-human CD11b, Ultra-LEAF™ Purified anti-human CD11b, TotalSeq™-D0161 anti-human CD11b, GMP PE/Cyanine7 anti-human CD11b, GMP PE anti-human CD11b, Spark UV™ 387 anti-human CD11b, GMP PerCP/Cyanine5.5 anti-human CD11b, GMP APC anti-human CD11b, GMP FITC anti-human CD11b, GMP APC/Fire™ 750 anti-human CD11b, Spark Blue™ 515 anti-human CD11b, Spark Violet™ 500 anti-human CD11b, Spark Red™ 718 anti-human CD11b (Flexi-Fluor™), GMP Pacific Blue™ anti-human CD11b, APC/Fire™ 810 anti-human CD11b, Brilliant Violet 750™ anti-human CD11b Antibody

## Product Data

---



Human peripheral blood leukocytes were stained with anti-human CD11b (clone ICRF44) APC/Fire™ 810 (left) or mouse IgG1, κ APC/Fire™ 810 isotype control (right).

For Research Use Only. Not for diagnostic or therapeutic use.

This product is supplied subject to the terms and conditions, including the limited license, located at [www.biolegend.com/terms](http://www.biolegend.com/terms) ("Terms") and may be used only as provided in the Terms. Without limiting the foregoing, BioLegend products may not be used for any Commercial Purpose as defined in the Terms, resold in any form, used in manufacturing, or reverse engineered, sequenced, or otherwise studied or used to learn its design or composition without express written approval of BioLegend. Regardless of the information given in this document, user is solely responsible for determining any license requirements necessary for user's intended use and assumes all risk and liability arising from use of the product. BioLegend is not responsible for patent infringement or any other risks or liabilities whatsoever resulting from the use of its products.

BioLegend, the BioLegend logo, and all other trademarks are property of BioLegend, Inc. or their respective owners, and all rights are reserved.

8999 BioLegend Way, San Diego, CA 92121 [www.biolegend.com](http://www.biolegend.com)  
Toll-Free Phone: 1-877-Bio-Legend (246-5343) Phone: (858) 768-5800 Fax: (877) 455-9587