

Spark PLUS B488™ Mouse IgG1, κ Isotype Ctrl Antibody

Catalog# / Size	402454 / 25 tests 402455 / 100 tests
Clone	MOPC-21
Regulatory Status	RUO
Isotype	Mouse IgG1, κ
Description	The MOPC-21 immunoglobulin has unknown specificity. The isotype of this antibody is mouse IgG1, κ. This antibody was chosen as an isotype control after screening on a variety of resting, activated, live, and fixed mouse, rat and human tissues.

Product Details

Antibody Type	Monoclonal
Host Species	Mouse
Formulation	Phosphate-buffered solution, pH 7.2, containing 0.09% sodium azide and BSA (origin USA)
Preparation	The antibody was purified by affinity chromatography and conjugated with Spark PLUS B488™ under optimal conditions.
Concentration	Lot-specific (to obtain lot-specific concentration and expiration, please enter the lot number in our Certificate of Analysis online tool.)
Storage & Handling	The antibody solution should be stored undiluted between 2°C and 8°C, and protected from prolonged exposure to light. Do not freeze.
Application	FC - Quality tested
Recommended Usage	Each lot of this antibody is quality control tested by immunofluorescent staining with flow cytometric analysis as negative control. Use at concentrations comparable to those of the specific antibody of interest. Use our Certificate of Analysis online tool to find the exact concentrations of your lots of product. * Spark PLUS B488™ has a maximum excitation of 499 nm and a maximum emission of 520 nm.
Excitation Laser	Blue Laser (488 nm)
Application Notes	Additional reported applications (for the relevant formats) include: Intracellular Flow Cytometry (ICFC), Immunocytochemistry (ICC), Immunohistochemistry (IHC), Immunoprecipitation (IP), Western Blotting (WB), Functional Assay (FA)
Application References (PubMed link indicates BioLegend citation)	1. Carlsten M, <i>et al.</i> 2007. <i>Cancer Res.</i> 67:1317. PubMed 2. Smed-Sørensen A, <i>et al.</i> 2008. <i>Blood</i> 111:5037. PubMed (FA) 3. Bunesmann MM, <i>et al.</i> 2011. <i>Am. J. Respir. Cell. Mol. Biol.</i> Epub. PubMed 4. Matsuyama T, <i>et al.</i> 2005. <i>Infect. Immun.</i> 73:1044. (IF) 5. Correia DV, <i>et al.</i> 2011. <i>Blood</i> 118:992. (FC) PubMed 6. Lian IA, <i>et al.</i> 2011. <i>Placenta.</i> 32:823. PubMed 7. Bufe B, <i>et al.</i> 20015. <i>J Biol Chem.</i> 290:7369. PubMed

Antigen Details

Related Protocols

- [Cell Surface Flow Cytometry Staining Protocol](#)

Other Formats

Spark Violet™ 538 Mouse IgG1, κ isotype Ctrl, Spark YG™ 581 Mouse IgG1, κ Isotype Ctrl, Spark YG™ 593 Mouse IgG1, κ Isotype Ctrl, GMP FITC Mouse IgG1, κ Isotype Ctrl, GMP PE/Dazzle™ 594 Mouse IgG1, κ Isotype Ctrl, GMP PE Mouse IgG1, κ Isotype Ctrl, GMP APC Mouse IgG1, κ Isotype Ctrl, GMP PerCP/Cyanine5.5 Mouse IgG1, κ Isotype Ctrl, PE/Fire™ 810 Mouse IgG1, κ Isotype Ctrl, Spark YG™ 570 Mouse IgG1, κ Isotype Ctrl, TotalSeq™-Bn0090 Mouse IgG1, κ Isotype Ctrl, Spark Blue™ 574 Mouse IgG1, κ Isotype Ctrl, PerCP/Fire™ 806 Mouse IgG1, κ Isotype Ctrl, Alexa Fluor® 660 Mouse IgG1, κ Isotype Ctrl, Spark PLUS UV395™ Mouse IgG1, κ Isotype Ctrl, PE/Fire™ 744 Mouse IgG1, κ Isotype Ctrl, Alexa Fluor® 594 Mouse IgG1, κ Isotype Ctrl, Spark UV™ 387 Mouse IgG1, κ Isotype Ctrl, Spark Red™ 718 Mouse IgG1, κ Isotype Ctrl, Spark Violet™ 423 Mouse IgG1, κ Isotype Ctrl, GMP Pacific Blue™ Mouse IgG1, κ Isotype Ctrl, GMP PE/Cyanine7 Mouse IgG1, κ Isotype Ctrl, GMP APC/Fire™ 750 Mouse IgG1, κ Isotype Ctrl, Alexa Fluor® 700 Mouse IgG1, κ Isotype Ctrl (ICFC), Purified Mouse IgG1, κ Isotype Ctrl, GolnVivo™ Purified Mouse IgG1, κ Isotype Ctrl, PE/Cyanine5 Mouse IgG1, κ Isotype Ctrl, FITC Mouse IgG1, κ Isotype Ctrl, APC/Fire™ 750 Mouse IgG1, κ Isotype Ctrl, APC Mouse IgG1, κ Isotype Ctrl, Biotin Mouse IgG1, κ Isotype Ctrl, PE/Cyanine7 Mouse IgG1, κ Isotype Ctrl, PE Mouse IgG1, κ Isotype Ctrl, PE/Dazzle™ 594 Mouse IgG1, κ Isotype Ctrl, Brilliant Violet 570™ Mouse IgG1, κ Isotype Ctrl, Brilliant Violet 605™ Mouse IgG1, κ Isotype Ctrl, Ultra-LEAF™ Purified Mouse IgG1, κ Isotype Ctrl, Brilliant Violet 510™ Mouse IgG1, κ Isotype Ctrl, Brilliant Violet 421™ Mouse IgG1, κ Isotype Ctrl, Spark Violet™ 500 Mouse IgG1, κ Isotype Ctrl, Brilliant Violet 650™ Mouse IgG1, κ Isotype Ctrl, Spark Blue™ 515 Mouse IgG1, κ Isotype Ctrl, PerCP/Fire™ 780 Mouse IgG1, κ Isotype Ctrl, Alexa Fluor® 647 Mouse IgG1, κ Isotype Ctrl (FC), TotalSeq™-C1471 Mouse IgG1, κ Isotype Ctrl Antibody, Alexa Fluor® 555 Mouse IgG1, κ Isotype Ctrl, Spark PLUS B488™ Mouse IgG1, κ Isotype Ctrl, APC Mouse IgG1, κ Isotype Ctrl (ICFC), APC Mouse IgG1, κ Isotype Ctrl (FC), PE Mouse IgG1, κ Isotype Ctrl (FC), FITC Mouse IgG1, κ Isotype Ctrl (FC), Alexa Fluor® 488 Mouse IgG1, κ Isotype Ctrl (ICFC), FITC Mouse IgG1, κ Isotype Ctrl (ICFC), PE Mouse IgG1, κ Isotype Ctrl, PerCP Mouse IgG1, κ Isotype Ctrl, PerCP/Cyanine5.5 Mouse IgG1, κ Isotype Ctrl, Pacific Blue™ Mouse IgG1, κ Isotype Ctrl, APC/Cyanine7 Mouse IgG1, κ Isotype Ctrl, APC/Fire™ 810 Mouse IgG1, κ Isotype Ctrl, PE/Fire™ 640 Mouse IgG1, κ Isotype Ctrl, Brilliant Violet 711™ Mouse IgG1, κ Isotype Ctrl, Brilliant Violet 785™ Mouse IgG1, κ Isotype Ctrl, TotalSeq™-B1471 Mouse IgG1, κ Isotype Ctrl, Spark PLUS B550™ Mouse IgG1, κ Isotype Ctrl, Alexa Fluor® 488 Mouse IgG1, κ Isotype Ctrl (FC), Alexa Fluor® 647 Mouse IgG1, κ Isotype Ctrl (ICFC), PE Mouse IgG1, κ Isotype Ctrl (ICFC), TotalSeq™-A0090 Mouse IgG1, κ isotype Ctrl, TotalSeq™-C0090 Mouse IgG1, κ isotype Ctrl, PE/Fire™ 700 Mouse IgG1, κ isotype Ctrl, Alexa Fluor® 700 Mouse IgG1, κ Isotype Ctrl, Brilliant Violet 750™ Mouse IgG1, κ isotype Ctrl, TotalSeq™-B0090 Mouse IgG1, κ isotype Ctrl, KIRAVIA Blue 520™ Mouse IgG1, κ Isotype Ctrl, Spark NIR™ 685 Mouse IgG1, κ Isotype Ctrl, Spark Blue™ 550 Mouse IgG1, κ Isotype Ctrl

For Research Use Only. Not for diagnostic or therapeutic use.

This product is supplied subject to the terms and conditions, including the limited license, located at www.biolegend.com/terms ("Terms") and may be used only as provided in the Terms. Without limiting the foregoing, BioLegend products may not be used for any Commercial Purpose as defined in the Terms, resold in any form, used in manufacturing, or reverse engineered, sequenced, or otherwise studied or used to learn its design or composition without express written approval of BioLegend. Regardless of the information given in this document, user is solely responsible for determining any license requirements necessary for user's intended use and assumes all risk and liability arising from use of the product. BioLegend is not responsible for patent infringement or any other risks or liabilities whatsoever resulting from the use of its products.

BioLegend, the BioLegend logo, and all other trademarks are property of BioLegend, Inc. or their respective owners, and all rights are reserved.

8999 BioLegend Way, San Diego, CA 92121 www.biolegend.com
Toll-Free Phone: 1-877-Bio-Legend (246-5343) Phone: (858) 768-5800 Fax: (877) 455-9587