

## Spark PLUS UV395™ anti-human FOXP3 Antibody

<b>Catalog# / Size</b>	320137 / 25 tests 320138 / 100 tests
<b>Clone</b>	206D
<b>Regulatory Status</b>	RUO
<b>Other Names</b>	Forkhead box protein P3, Scurfin, JM2, IPEX, Zinc finger protein JM2
<b>Isotype</b>	Mouse IgG1, κ
<b>Description</b>	<p>FOXP3 is a 50-55 kD transcription factor, also known as Forkhead box protein P3, Scurfin, JM2, or IPEX. It is proposed to be a master regulatory gene and more specific marker of T regulatory cells than most cell surface markers (such as CD4 and CD25). Transduced expression of FOXP3 in CD4<sup>+</sup>/CD25<sup>-</sup> cells has been shown to induce GITR, CD103, and CTLA4 and impart a T regulatory cell phenotype. FOXP3 is mutated in X-linked autoimmunity-allergic dysregulation syndrome (XLAAD or IPEX) in humans and in "scurfy" mice. Overexpression of FOXP3 has been shown to lead to a hypoactive immune state suggesting that this transcriptional factor is a central regulator of T cell activity. In human, unlike in mouse, two isoforms of FOXP3 have been reported: one (FOXP3) corresponding to the canonical full-length sequence; the other (FOXP3 δ2) lacking exon 2. The 206D antibody recognizes human FOXP3 epitope in the region of amino acids 105-235.</p>

### Product Details

---

<b>Verified Reactivity</b>	Human
<b>Reported Reactivity</b>	Baboon, Cynomolgus, Rhesus, Pigtailed Macaque
<b>Antibody Type</b>	Monoclonal
<b>Host Species</b>	Mouse
<b>Immunogen</b>	Full-length FOXP3 protein
<b>Formulation</b>	Phosphate-buffered solution, pH 7.2, containing 0.09% sodium azide and BSA (origin USA)
<b>Preparation</b>	The antibody was purified by affinity chromatography and conjugated with Spark PLUS UV395™ under optimal conditions.
<b>Concentration</b>	Lot-specific (to obtain lot-specific concentration and expiration, please enter the lot number in our <a href="#">Certificate of Analysis</a> online tool.)
<b>Storage &amp; Handling</b>	The antibody solution should be stored undiluted between 2°C and 8°C, and protected from prolonged exposure to light. <b>Do not freeze.</b>
<b>Application</b>	<a href="#">ICFC – Quality tested</a>
<b>Recommended Usage</b>	<p>Each lot of this antibody is quality control tested by intracellular flow cytometry using our True-Nuclear™ Transcription Factor Staining Protocol. For flow cytometric staining, the suggested use of this reagent is 5 μL per 10<sup>6</sup> cells in 100 μL volume. It is highly recommended that the reagent be titrated for optimal performance for each application.</p> <p>* Spark PLUS UV395™ has a maximum excitation of 355 nm and a maximum emission of 385 nm.</p>
<b>Excitation Laser</b>	Ultraviolet Laser (355 nm)
<b>Application Notes</b>	Additional reported applications (for the relevant formats) include: immunohistochemical staining of acetone-fixed frozen sections <sup>1</sup> and formalin-fixed paraffin-embedded sections <sup>1,8,19-20</sup> , and Western blotting <sup>1</sup> . The binding of 206D to FOXP3 can be partially blocked by 259D, but 206D does not show significant blocking effect on 259D binding.

**NOTE:** For flow cytometric staining with this clone, True-Nuclear™ Transcription Factor Buffer Set (Cat. No. [424401](#)) offers improved staining and is highly recommended.

#### Application References

(PubMed link indicates BioLegend citation)

1. Roncador G, *et al.* 2005. *Eur. J. Immunol.* 35:1681.(IHC)
2. Yang ZZ, *et al.* 2006. *Blood* 107:3639.
3. Liu W, *et al.* 2006. *J. Exp. Med.* 203:1701.[PubMed](#)
4. Bollyky PL, *et al.* 2007. *J. Immunol.* 179:744.
5. Bell MP, *et al.* 2007. *J. Immunol.* 179:1893.
6. Tran DQ, *et al.* 2007. *Blood* doi:10.1182/blood-2007-06-094656. [PubMed](#)
7. Gao Q,*et al.*2007.*J Clin Oncol.*25:2586.(IHC) [PubMed](#)
8. Pillai V,*et al.* 2008. *Blood* 111:463. [PubMed](#)
9. Zheng Y, *et al.* 2008. *J. Immunol.* 181:1683. [PubMed](#)
10. Zonios DI, *et al.* 2008.*Blood*112:287. [PubMed](#)
11. Kavanagh B, *et al.* 2008. *Blood.* [PubMed](#)
12. Nevala WK, *et al.* 2009. *Clin Cancer Res.* 15:1931. [PubMed](#)

[See More](#)

#### RRID

AB\_3683337 (BioLegend Cat. No. 320137)  
AB\_3683337 (BioLegend Cat. No. 320138)

## Antigen Details

---

<b>Structure</b>	Forkhead/winged-helix transcription factor family, approximately 50 kD, contains zinc finger and forkhead domains
<b>Distribution</b>	Nuclear; expressed in T regulatory cells
<b>Function</b>	Transcription factor proposed to be a master regulatory gene in T regulatory cell development and a critical factor for immune homeostasis
<b>Interaction</b>	Interacts with DNA
<b>Cell Type</b>	Tregs
<b>Biology Area</b>	Cell Biology, Immunology, Transcription Factors
<b>Molecular Family</b>	Nuclear Markers
<b>Antigen References</b>	1. Hori S, <i>et al.</i> 2003. <i>Science</i> 299:1057. 2. Gandhi R, <i>et al.</i> 2010. <i>Nat. Immunol.</i> 11:846.
<b>Regulation</b>	FOXP3 is present at high levels in T regulatory cells, it can also be induced by T cell activation.
<b>Gene ID</b>	<a href="#">50943</a>

## Related Protocols

---

- [Intracellular Flow Cytometry Staining Protocol](#)


## Other Formats

---

Purified anti-human FOXP3, Alexa Fluor® 488 anti-human FOXP3, Alexa Fluor® 647 anti-human FOXP3, FITC anti-human FOXP3, Pacific Blue™ anti-human FOXP3, PE anti-human FOXP3, PE/Dazzle™ 594 anti-human FOXP3, True-Nuclear™ One Step Human Treg Flow Kit, Brilliant Violet 421™ anti-human FOXP3, KIRAVIA Blue 520™ anti-human FOXP3, Spark NIR™ 685 anti-human FOXP3 Antibody, Brilliant Violet 650™ anti-human FOXP3, Brilliant Violet 785™ anti-human FOXP3, Spark PLUS UV395™ anti-human FOXP3

## Product Data

---



Human peripheral blood lymphocytes were treated with True-Nuclear™ Transcription Factor Buffer Set and were then surfaced stained with anti-human CD4 (clone RPA-T4) APC and intracellularly stained with anti-human FOXP3 (clone 206D) Spark PLUS UV395™ (left) or mouse IgG1, κ Spark PLUS UV395™ isotype control (right).

For Research Use Only. Not for diagnostic or therapeutic use.

This product is supplied subject to the terms and conditions, including the limited license, located at [www.biolegend.com/terms](http://www.biolegend.com/terms) ("Terms") and may be used only as provided in the Terms. Without limiting the foregoing, BioLegend products may not be used for any Commercial Purpose as defined in the Terms, resold in any form, used in manufacturing, or reverse engineered, sequenced, or otherwise studied or used to learn its design or composition without express written approval of BioLegend. Regardless of the information given in this document, user is solely responsible for determining any license requirements necessary for user's intended use and assumes all risk and liability arising from use of the product. BioLegend is not responsible for patent infringement or any other risks or liabilities whatsoever resulting from the use of its products.

BioLegend, the BioLegend logo, and all other trademarks are property of BioLegend, Inc. or their respective owners, and all rights are reserved.

8999 BioLegend Way, San Diego, CA 92121 [www.biolegend.com](http://www.biolegend.com)  
 Toll-Free Phone: 1-877-Bio-Legend (246-5343) Phone: (858) 768-5800 Fax: (877) 455-9587